



Rev 03, EN, Quick Refer

9952.



Sling Positioning

Applying sling from a seated position




Stand behind or to the side of the patient and hold the sling with the label facing up and away from the patient. Have patient lean forward (assist if needed).



Place sling behind patient's back. The bottom of the sling should be in line with patient's coccyx. Pull the leg sections around the patient's hips.





their skin.

leg.

Roll the patient on their back and gently pull the gathered towards the bed to prevent skin

Pull the leg straps around the

Stand to side of patient and turn

them to their side. Place the

gathered sling with the label

and facing up and away from

towards the patient's head end



Positioning the legs



Open position: To achieve this This position is not available for this type of sling

Divided position: To achieve

This position is not

rap

r in

he

aps

rts

nd

available for this

type of sling

thi

su

thi


thi

аc



Positioning the body





Important! Shoulder loops should always be attached before the leg loops. Attaching the leg loops first may cause risk of the patient sliding out if the lifting process is started and all the straps are not secure.

Upright position: To achieve

this position, use the short

leg straps.

loops of the head and shoulder

straps and the long loops of the

Reclined position: To achieve

this position, use the middle

or long loops of the head and

shoulder straps and the middle

or short loops of the leg straps.

🗑 Important! The leg straps should surround the hips and be pulled to the same length on both sides! Colour match the loops! Do not let the material fold under the patient!




Lower the lift to the height of the patient's chest and affix the upper straps followed by the leg straps, ensuring colours are matched. Hold the lift with one hand away from the patient and slowly raise.



Applying sling from a supine position

Important! The leg straps should surround the hips and be pulled to the same length on both sides! Colour match the loops! Do not let the material fold under the patient!



Closed position: To achieve this position, pull one leg support, then the other under both thighs and attach to the end of the carry bar.

Order of strap attachment



Models: 25064, 25069 Human Care HC AB. All Rights Reserved.





Lower the lift to the height of the patient's chest and affix the upper straps followed by the leg straps, ensuring colours are matched. Hold the lift with one hand away from the patient and slowly raise.