



Table Pad  
Indoor / Outdoor  
Unique Handles



The Rebel Trolley is a safe and reliable rollator with a firm table pad for extra support. Designed for users requiring that extra support, or those who find it difficult to maneuver an ordinary walker, the Rebel TW is the ideal rollator for both indoor or outdoor use.

#### Outstanding Stability

The Rebel TW is the ideal rollator for pain sensitive users suffering from arthritis, Parkinson's or many other neurological disorders. It features a wide, firm table support and large ball-bearing mounted wheels that ensure easy maneuverability and comfort. The Rebel TW can be used both indoors and outdoors and features a wide assortment of accessories for everyday use.

#### Unique Handles

The Rebel TW features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle

tube. The ergonomic design of the grips ensures even distribution of weight. The Rebel TW can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

#### Reactive Frame

The Rebel TW has the same great frame as the regular Rebel which means that it remains stable even on uneven surfaces. This is especially helpful in rural areas where the ground is not even. The Rebel TW comes in two seat heights to accommodate most users.

#### HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B  
117 43 Stockholm  
Phone: +46 8 665 35 00  
Fax: +46 8 665 35 10  
info.se@humancaregroup.com

#### HUMAN CARE CANADA

10-155 Colonnade Road  
Ottawa, ON K2E 7K1  
Phone: +613.723.6734  
Fax: +613.723.1058  
info.ca@humancaregroup.com

#### HUMAN CARE NEDERLAND

Elspeterweg 124  
8076 PA Vierhouten  
Phone: +31 577 412 171  
Fax: +31 577 412 170  
info.nl@humancaregroup.com

#### HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K  
Austin, TX 78754  
Phone: +512.476.7199  
Fax: +512.476.7190  
info.us@humancaregroup.com

## Accessories/Features



### Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10601



### Basket

Included with rollator in Canada. Ideal for shopping. Part no. 10621



### Reflector Kit

Includes 30 self-adhesive reflectors. Part no. 91090



### Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620



### Light kit

Illuminate your way while walking. Part no. 91091



### Security Lock

Flex cord that protects rollator from theft while unattended. Part no. 91092




### Cane holder

Simplifies transport of a cane or umbrella. Part no. 10603



### One-Hand Brake

Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10653



### Oxygen Holder

Holds oxygen tube firmly in place. Part no. 10615



### Curb Climber

Helps raise rollator onto a curb. Part no. 10612

Note: Not all accessories are available for sale in all countries.

## Technical Specifications

	Rebel TW Low 67	Rebel TW Low 72	Rebel TW Standard 67	Rebel TW Standard 72
Part no.	23311B	23411B	23321B	23421B
Handle Height (Min - Max)	86cm - 104cm / 33.8" - 40.9"	86cm - 104cm / 33.8" - 40.9"	104cm - 122cm / 40.9" - 48"	104cm - 122cm / 40.9" - 48"
Width between Handles	40.6cm / 16"	40.6cm / 16"	50.8cm / 20"	50.8cm / 20"
Wheel Size	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"
Product Weight	12.8kg / 28.2lb	12.9kg / 28.4lb	13.0kg / 28.7lb	13.1kg / 28.9lb
Maximum User Weight	125kg / 275.6lb	125kg / 275.6lb	125kg / 275.6lb	125kg / 275.6lb