



Easy to fold  
Reactive frame  
Unique handles



The Carl-Oskar is the perfect rollator for the active user. Its durable frame is easy to maneuver and features ergonomic handles for a comfortable grip. With extra large wheels, it is the perfect rollator for those outdoor activities, even during snowy winters.



Carl-Oskar  
Reinforced

### Outstanding Stability

The Carl-Oskar is one of Human Care's best known and best-selling rollators worldwide. It is safe and reliable and features extra-large 30cm/11" wheels that allow the rollator to roll very easily yet remain stable on virtually any surface. The Carl-Oskar is the ideal option for rural areas, even on snowy surfaces.

functions on either the left or right hand side.

### Reinforced Option

The Carl-Oskar is available in different overall widths to accommodate a wide variety of users. The reinforced models for bariatric users feature a weight capacity of 275kg/606lb. All models come with a full assortment of accessories to meet most needs and can easily be folded for storage or transport.



### Unique Handles

The Carl-Oskar features unique angle adjustable handles that allow users to find the optimal handle position for comfort and stability. The ergonomic design ensures even distribution of weight for a stable grip. The Carl-Oskar can also be equipped with a one-hand brake to allow a user to control hand

#### HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B  
117 43 Stockholm  
Phone:- +46 8 665 35 00  
Fax:- +46 8 665 35 10  
info.se@humancaregroup.com

#### HUMAN CARE CANADA

10-155 Colonnade Road  
Ottawa, ON K2E 7K1  
Phone: +613.723.6734  
Fax: +613.723.1058  
info.ca@humancaregroup.com

#### HUMAN CARE NEDERLAND

Elspeterweg 124  
8076 PA Vierhouten  
Phone: +31 577 412 171  
Fax: +31 577 412 170  
info.nl@humancaregroup.com

#### HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K  
Austin, TX 78754  
Phone: +512.476.7199  
Fax: +512.476.7190  
info.us@humancaregroup.com

## Accessories/Features



**Tray**  
Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10651



**Splash Guard**  
Protects legs and clothes. Part no. 10608



**Cross handle**  
Cross handle with included one-hand brake. Only for Carl-Oskar 62. Art.no: 70844



**Slow-Down Brake**  
Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10624



**Light kit**  
Illuminate your way while walking. Part no. 91091



**Security Lock**  
Flex cord that protects rollator from theft while unattended. Part no. 91092



**Cane holder**  
Simplifies transport of a cane or umbrella. Part no. 10603



**One-Hand Brake**  
Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10652



**Bag Kit**  
Includes one small (5kg capacity) and on large bag (10kg capacity). Part no. 91095 (for 62cm), 91096



**Oxygen Holder**  
Holds oxygen tube firmly in place. Part no. 10615



**Seat Pad**  
State-of-the-art technology keeps pad cool in summer and warm in winter. Part no. 91093



**Backstrap and Basket**  
Included with rollator in Canada. Backrest Part no. 10621, 10655 (for 72cm) Basket Part no. 10605

## Technical Specifications

Description	Part No	Width between Handles	Overall Width	Distance btwn Rear Wheels	Product Weight
Carl-Oskar 57	11150 11151 (Air Tires)	37cm / 14.6"	57cm / 22.4"	50.8cm / 20"	11.1kg / 24.5lb
Carl-Oskar 62	11250 11251 (Air Tires) 61225 (Reinforced)	42cm / 16.5"	62cm / 24.4"	50.8cm / 20"	11.3kg / 24.9lb 11.3kg / 24.9lb 13.2kg / 29.1lb
Carl-Oskar 67	11350 11351 (Air Tires) 61325 (Reinforced)	47cm / 18.5"	67cm / 26.4"	52.1cm / 20.5"	11.4kg / 25.1lb 11.4kg / 25.1lb 13.3kg / 29.3lb
Carl-Oskar 72	11450 11451 (Air Tires) 61425 (Reinforced)	52cm / 20.5"	72cm / 28.3"	52.1cm / 20.5"	11.5kg / 25.3lb 11.5kg / 25.3lb 13.4kg / 29.5lb

Note: Not all accessories are available for sale in all countries.

Description (All Models)	
Handle Height (Min - Max)	77cm - 95cm / 30.3" - 37.4"
Seat Height	64cm / 25.2"
Wheel Size	30cm x 5cm / 11.8" x 1.8"
Seat Size	31.8cm x 17.8cm / 12.5" x 7"
Depth folded	26cm / 10.2"
Maximum User Weight	Regular 150kg / 330.7lb Reinforced 275kg / 606.3lb