

# Multi Sling Fitting Protocol



**Client Information:** Name:  Room:

Facility:  Lift Used:  Hanger Bar Used:

**Sling Model:** Name/Article Number:  Fabric:  Size:

**Loops for back:**  **Green Strap**  **Yellow Strap**  **Red Straps**

**Extension Loops Back:**

21180-10 [10cm/3.9"]
21180-15 [15cm/5.9"]
21180-20 [20cm/7.8"]
21180-25 [25cm/9.8"]
21180-35 [35cm/13.7"]

\*Please note that the leg placement information box below is not relevant for the Multi Sling. \*Only arm placement on the outside of the Sling possible for the Multi Sling. \*Only 2-point and 4-point suspension is possible for the Multi Sling.

**Transfer Situation:**

**Arm placement:**

**Leg Placement:**



**Suspension:**

**Accessories required:**  50866 Upholstered Cover for Hanger Bar   
**Do not use in wet environments.**

**Notes:**

Fitting Protocol Completed by:  Position/Title:  Date:

**Size guidelines:**

 User weight kg	 User length cm	 Waist circumfer- ence cm	 Recommended size*	 Color
30-60	140-170	60-91	Small (S)	
55-90	160-190	94-110	Medium (M)	
85-120	180-206	104-129	Large (L)	
110-200	180-210	129-145	Extra Large (XL)	

\* These sizes are recommendations only. A correct fitting must always be carried out by clinically knowledgeable personnel using the sling's user manual. Please be aware that some slings are not available in all sizes. You can find drawings on each sling in the user manual, delivered with the product and available on our websites.

Fitting Protocol Completed by:

Name:

Date:  Position/Title:



**Important!** This fitting protocol does not replace the user manual. Read all instructions before use!