


Alpha Soft Turning Disc

Simple
Soft
Sideward
rotation



The Alpha Soft Turning Disc is a simple aid that significantly facilitates lateral movements for seated clients, whether alone or with trusted help.

For use in a bed, place the turning disc centrally in the bed approx. 25 cm from the edge of the bed. You thus achieve an optimal lying position in the centre of the bed from the turning movement.

While in the seating position on the bed always use the height adjustment of the bed to assist with your safety.

Our Soft Turning Aid is made of breathable cotton.

Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Soft Turning Disc	PM-5070	46 cm diameter / 18,1" diameter	Cotton/latex free	200 kg / 440 lb



PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.