



Non-satic
Head support
Grips



The Alpha Caremaster sliding mat is used for the transfer of clients in seated and supine positions and for turning patients in bed. It assists by providing secure support and momentum to the movement sequence without having to touch the client in a sensitive area.

The mat has a soft padding to allow comfortable transfer and the integrated head support offers the possibility to perform a gentle transfer. The three extra strong webbed grips can also help when handling very heavy persons.



Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Caremaster Medium	PM-2095	60 x 175 cm 68,9" x 23,6"	PES/Carbon, latex free	200 kg 440 lb
Alpha Caremaster Extra Large	PM-2097	70 x 185 cm 72,8" x 27,6"	PES/Carbon, latex free	200 kg 440 lb
Alpha Caremaster Extra Extra Large	PM-2094	220 x 80 cm / 86,6 x 31,5"	PES/Carbon, latex free	200 kg 440 lb



PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.