


Alpha Flexi Board Oval

Smooth
Easy to use
Flexible



Ergonomically designed handles.



Easy to hang and store.

Alpha Flexi Board Oval is available in two different sizes and is used to move and reposition users.

Alpha Flexi Board Oval L is used for horizontal transfers, positioning in bed and application of slings. The sliding board has low friction on both sides and evens out unstable surfaces which benefits a soft surface, and facilitates movements.

Alpha Flexi Board Oval S is used to get back in a (wheel)chair, minor repositioning and when applying slings.

The oval shape of the sliding board together with the thin low-friction material makes it possible to easily and smoothly place and remove the board from the user.

The handles are ergonomically designed and allow the caregiver to get a good grip during use and movement. The holes on each long side are intended as an alternative hanging option for the sliding board, eg near the bed, so the product is always ready for use.

Product specification

Product name	Art. nr.	Dimensions	Material
Alpha Flexi Board Oval L (2 pieces)	PM-5040	85 x 40 cm 33,5" x 15,7"	Polyethene, latex free
Alpha Flexi Board Oval S (2 pieces)	PM-5045	46 x 20 cm 18,1" x 7,9"	Polyethene, latex free



 = 150 kg / 330 lb

