

Guess 402

Ergonomic
Adjustable
Safe



The comfortable, cost-effective Guess 402 is suitable for daily patient care, as well as for easier intensive care and to promote rapid recovery. Its modern control options are particularly attractive. Guess 402 split gates are easy to use and have excellent hygienic and functional properties.

The hospital bed that meets all requirements

Hospital beds are expected to meet an extremely wide range of requirements. The tasks they are faced with at the various hospital units are often extensive, which means that several different models must be purchased.

We know that the market has a need for beds with many different functions adapted to the unit where the bed is to be used. Seen from a user, environmental and cost perspective, it would be more efficient to use a bed model that meets the requirements that different care units have for their operations. This is exactly the approach that the GUESS 402 from KENMAK follows.

Its smart modular system makes the GUESS402 competent in all hospital units. It is equipped with four motors, and a variety of control options, safety systems and accessories which make it possible to create a practical, attractive and versatile hospital bed that meets individual requirements while being cost-effective.

Guess 402 makes the work of the staff much easier by being ergonomically designed and easy to use. The design of the bed also facilitates cleaning, which is important from the hygiene requirements set within the operations.

Kenmak
part of the Human Care Group

Human Care AB is a distributor
of Kenmak



Walking Aids



Lifting Solutions



Healthcare Beds



Hygiene



Convertible chairs

Accessories



Integrated scale

Sensors sit in each corner of the bed base and measure weight automatically, 0-400kg.

Art nr: G402SCALE



Lifting Pole: Art nr: G402LP

IV-Pole: Art nr: G402LVP

IVA-Pole available Art nr: G402IVA



Monitor Shelf

Stable table as folding. Placed on head or footboard.

Art nr: G402MS

Document holder

Art nr: G402DH



Integrated siderail control panels

For this bed we also offer integrated caregiver control in the bed gate. Art nr: G402IP



5th wheel

Facilitates maneuvering of the bed in tight spaces.

Art nr: G4025THWCOMPL



Nurse control panel

Battery indicator, locking functions for buttons. Nursing staff has one-button functions, button for CPR and Shock mode.



Remote control panel

Control for both patient and healthcare staff.



Trendelenburg

16 degrees Trendelenburg / reverse Trendelenburg



Head and Foot Boards

Easily removed without tools. Made of ABS plastic.



Side Rails Gates with support handles in ABS plastic. The bed is equipped with two gates on each side.

Siderail extension is available



Bed Extension

Integrated extension of bed: 20 cm.

Art nr: G402BEEXT



Central locking

Antistatic, 150 mm encapsulated double wheels with central locking. More sizes of wheels are available, see our website.

Technical specifications

Guess 402

Art. nr:	Guess 402-SE	Height Range (Min-Max)	38 cm - 80 cm / 15 - 31,5" (without mattress)
Sleep Plane (L x W)	200 cm x 85 cm / 78,8 x 33,5"	Bed Measurements (L x W)	220 cm x 104 cm / 86,6 x 50"
Rated Voltage Distance	220 / 240 Volt	IP-number	IPX6
Rated Frequency	50/60Hz	Product Weight	120 kg / 264.6 lb
Maximum User Weight	185 kg / 408 lb	Max load capacity	250 kg / 551 lb
Mode of operation	2 minutes on / 18 minutes off		

